




The design features:

- **COMFORTABLE** three finger ring grip for ease of injection of viscid Phenol.
- **ANGLED or STRAIGHT** luer lock needle for security and improved ease of access.
- **WIDE BORE** needle with 18G distal tip to provide maximum control during injection.
- **SUPPLIED COMPLETE** with filling cannula for ease of use in the busy clinic setting.
- **RELIABLE**, practical and competitively priced to ensure value for your department.

Haemorrhoid Injection Sets provide a safe, reliable method of injection of Oily Phenol BP, or other sclerosing agents, as an outpatient treatment for haemorrhoids.



Ordering Information		NHS	Code
<b>Haemorrhoid Injection Set</b> Supplied with 10ml, 3 ring syringe, haemorrhoid injection needle, filling cannula. Supplied sterile, for single use. 25 sets per carton.	Straight, 10cm, 4"	FTK063	<b>R56210</b>
	Angled, 10cm, 4"	FTK011	<b>R56211</b>