

- **ROCKET MEDICAL SPINAL MANOMETERS WITH LUER FITTINGS** connect to luer lock spinal needles and syringes to facilitate the removal and measurement of Cerebrospinal Fluid (CSF).
- **READY FOR IMMEDIATE USE** for when time is precious.
- **CLEAR, STURDY PLASTIC CONSTRUCTION** for reliable, secure assembly and ease of use.
- **STACKABLE ASSEMBLY:** Male/female luer tube fittings permit additional sections to be stacked, allowing pressures >40cmH₂O to be measured when required.
- **SUPPLIED COMPLETE** with 3-way stopcock for attachment to all luer fitting spinal needles and syringes.



Rocket Medical Spinal Manometers have been providing safe, reliable, cost effective CSF pressure measurement for over 30 years.

The unique 'stackable' design, original to Rocket, has become a standard neurological diagnostic tool and can be found in emergency departments, neurological wards, medical wards and specialist neurosurgery centres throughout the world.

Description	Order Code	NHS SC Code
Spinal Manometer - Luer Lock 1 x Spinal Manometer Set with 3 way stopcock. Sterile, for single use, packed in cartons of 40 units	R55990	FTP002