



Rocket Medical lightweight aspiration pump filter sets offer a simple, safe, cost effective, patient connection method specifically designed for use with the latest generation of digital Oocyte Aspiration Pumps.

- ◆ **REMOVES CONTAMINATION RISK** from irregular or uncontrolled decontamination of multiple use consumable filters.
- ◆ **LIGHTWEIGHT PATIENT TUBE SET:** 2.5m lightweight patient tube set with integral luer lock connector to ensure secure, reliable connection to all standard luer-lock fitting aspiration needle sets.
- ◆ **HYDROPHOBIC BARRIER** provides pump protection by preventing accidental aspiration of fluid into the device.
- ◆ **DESIGNED SPECIFICALLY** to match the performance characteristics of digital Oocyte Aspiration Pumps

FILTER SPECIFICATION

- 17mm
- Effective filtration 0.2micron,
- BFE 99.9999%*
- VFE 99.9999%*
- High flow /low resistance Insufflation/Gas/Air Filter
- Hydrophobic barrier.
- Inner core: PTFE

Description	Order Code
Filter Set for Oocyte Aspiration Pump 0.2µm filter set, 2.5m patient connection tube with male luer connector. Supplied sterile, for single use in cartons of 10 units.	R57686
Filter Set for Cook Medical K-MAR-5200 Pumps 0.2µm filter set, 2.5m patient connection tube with male luer connector. Supplied sterile, for single use in cartons of 10 units.	R57687