



# INTRAUTERINE SOUND



- **WELL TOLERATED:** 3.5mm OD distal bead and 2.5mm distal shaft, for easier, low-trauma insertion.
- **MALLEABLE DISTAL SHAFT** allows for increased angulation, for easier insertion into the retroverted or anteverted uterus.
- **GRADUATED MARKINGS:** 1cm depth markers commence 40mm from the distal tip. The proximal marker is positioned 90mm from the distal tip.
- **HANDLE DESIGN** increases grip and tip control, whilst reducing bulk to improve visualisation of the cervix.
- **STERILE** and individually wrapped, for single patient use. Supplied in cartons of 100 units.

Description	Order Code
<b>Intrauterine Sound</b> Sterile, for single patient use. Packed in cartons of 100 units	<b>R57420</b>

