



- **READY TO USE:** self-contained battery powered light source with a ring of six high intensity white light LED's for maximum illumination.
- **MAGNIFYING LENS** for clear and precise vision and easy application of biopsy instrumentation
- **COMPLETE SET** supplied with a disposable insufflation bulb and pilot balloon to provide a complete and economical single use system.
- **SAFETY:** Should local hospital policy require it, the battery pack is fully removable to permit safe, environmentally responsible disposal. Each carton is supplied with a battery storage tray to ensure safe separation of used battery cells.

The Rocket® Pro-Lux™ fully disposable sigmoidoscopy system is a ready to use device equipped with new features to make usage convenient, practical and effective.



Description	Order Code	CH2 Code
Rocket® Pro-Lux™ Sigmoidoscope Self-contained illumination handle, insufflation bulb, pilot balloon. Non-sterile, 25 units per carton	R53901	2129924
Rocket® Pro-Lux™ Sigmoidoscope - Lubricated Self-contained illumination handle, insufflation bulb, pilot balloon. Non-sterile, 25 units per carton.	R53901-L	2129932


CLIFFORDHALLAMHEALTHCARE

 1300 720 274 www.ch2.net.au

nusurgix is a trading name of Rocket Medical plc

Rocket Medical plc
Suite 6,157 Gordon Street, Port Macquarie NSW 2444. AUSTRALIA.
Tel: +61 (0) 431 146448 Email: aus@rocketmedical.com W: aus.rocketmedical.com

 SALES 237 Rev: 03
 Date: 2021-06-15
 Copyright© 2014-21