

FBS Kit for Handheld Lactate Analysers



- ◆ **EASIER SAMPLE COLLECTION:** The sampling wand combines blade and capillary tube holders, to facilitate continuous visibility of the sampling site, which can be lost when transferring between separate devices.
- ◆ **SAFER COLLECTION:** The sampling wand features a retractable blade to reduce the potential for needlestick injuries. To retract the blade after use, just press the button on the shaft - it is as easy to use as a ballpoint pen.
- ◆ **MAXIMUM ILLUMINATION AS STANDARD:**
To ensure bright and reliable visualisation, every kit contains the latest **Amnilume™ light source**.
- ◆ **CAPILLARY TUBE DISPENSERS:** Each kit contains three individual capillary tube dispensers. Simply remove the filled capillary tube (distal 1-2cm) from the Sampling Wand or capillary tube holder, insert into the adaptor on the tube dispenser bulb and lightly depress the bulb to deliver the blood sample directly to the test strip.



At Rocket Medical we know our customers are our greatest source of new ideas and product improvements.

So, our latest FBS Kits for handheld lactate monitors incorporate 3 dedicated capillary tube dispensers to make delivering the low volume sample typically used by strip based monitors, quicker, easier and safer.

CONTENTS:

- Amnioscope and Amnilume™ Light Source:
142mm long x 33mm distal OD
- Sampling Wand - combined blade and capillary tube holder
- Double-ended capillary tube holder
- 5 x 10cm x-ray detectable gauze swabs - green.
- 6 x swab sticks
- 2.5gms petroleum jelly
- Sponge-holding forceps
- 5 x 55µl, high-heparin capillary tubes
- tube covers
- 3 x capillary tube dispensers
- Packed in cartons 10 units, for single use.

FBS Kit for handheld lactate monitors, 5-10µl sample volume

Code

FBS Kit for Handheld Lactate Analysers

R57028-00-SW

Packed in cartons 10 units, sterile, for single use