



Embryon[®] BulbTip[™] Embryo Transfer



Embryon[®] BulbTip[™] & Embryon[®] BulbTip[™] Echogenic catheters

- **ULTRAFINE INNER TRANSFER CATHETER** to minimise trauma and media volume on transfer. Externally supported to ease loading and to aid insertion into the sheath during transfer.
- **BULB TIP OUTER SHEATH** with smooth bulb formed tip to improve location during insertion, reduce cervical trauma and to guide the ultrafine inner catheter correctly into the uterine cavity. Particularly suitable for patients where cervical access is difficult.
- **SLIDING SILICONE MARKER** to allow insertion depth to be pre-set prior to transfer.
- **ECHO BANDED VERSION** to improve visualisation under ultrasound guidance.
- **RIBBED FINGER GRIP** gives unequalled grip and control while ensuring accurate placing of the catheter tip into the uterine cavity.



Embryon[®] BulbTip[™] & Embryon[®] BulbTip[™] Echogenic catheters are the latest generation of embryo transfer products from Rocket Medical.

Designed with innovative features to further improve ease of insertion and minimize uterine trauma, they continue our long established reputation for continuous product development.

Description	Length	Order Code
Embryon[®] BulbTip[™] Embryo Transfer Catheter Sterile, for single use, supplied in cartons of 10 units	18cm	R57635-00-18
	23cm	R57635-00-23
Embryon[®] BulbTip Echogenic Embryo Transfer Catheter Sterile, for single use, supplied in cartons of 10 units	18cm	R57635-EC-18
	23cm	R57635-EC-23
Embryon[®] BulbTip[™] Trial Embryo Transfer Catheter Sterile, for single use, supplied in cartons of 10 units	18cm	R57636-00-18
	23cm	R57636-00-23