




- **TWO PART CONSTRUCTION IUI** with smooth outer sheath with high memory characteristics ensures ease of insertion, particularly in those patients where access to the external os is difficult.
- **SOFT INNER IUI** with fully domed distal tip and dual side eyes for optimum insemination.
- **COMPETITIVELY PRICED** to provide excellent value in routine IUI.
- **STERILE FOR SINGLE USE**  
- Packed in cartons of 25 units.



**Rocket Medical DUO™ IUI Catheters** are the latest generation of two part insemination catheters designed to improve ease of insertion and minimise uterine trauma.

Used throughout the world, Rocket Medical DUO™ Catheters are the device of choice when access is restricted and insemination likely to be difficult.

**Specifications:**

- (a) Inner Transfer Catheter: Polyurethane OD: 1.4 mm, ID: 0.9 mm, fully domed, dual side eye.  
 (b) Outer Sheath: HDPE, 5 x 1cm marking from distal tip.  
 (c) Inner Catheter volume: 18cm: 0.5ml, 23cm: 0.63ml.

Description	Length	Order Code
<b>DUO™ IUI Catheters.</b> Sterile, for single use, packed in cartons 25 units	18cm	<b>R57625</b>
	23cm	<b>R57625-00-23</b>