



- **EDL™ DOUBLE LUMEN** catheter with integral stylet wire, 5 x 1cm marking, atraumatic distal tip and dual side eyes for insemination where a malleable catheter is required to gain access.
- **COMPETITIVELY PRICED** to provide excellent value in routine IUI.
- **STERILE FOR SINGLE USE** - Packed in cartons of 25 units.

**Rocket Medical EDL™ Double Lumen IUI Catheter**, includes an integral stiffening wire for insemination where cervical access is complicated by excessive retroversion, anteversion or cervical stenosis.



Tip detail of EDL™ IUI Catheter

**Specifications:**

**EDL™**

Outer: 133 x 2.60 x 2.23mm. - polypropylene

Inner: 173 x 2.10 x 0.55mm

40mm active length – polythene.

Stylet wire: stainless steel.

Description	Length	Order Code
<b>EDL™ IUI Catheter</b> Sterile, for single patient use, peel pouch, packed in cartons 25 units	18cm	<b>R57622</b>