

Rocket® Pro-Lux™ Sigmoidoscope



The Rocket® Pro-Lux™ fully disposable sigmoidoscopy system is a ready to use device equipped with new features to make usage convenient, practical and effective.

- **READY TO USE:** self-contained battery powered light source with a ring of six high intensity white light LED's for maximum illumination.
- **MAGNIFYING LENS** for clear and precise vision and easy application of biopsy instrumentation
- **COMPLETE SET** supplied with a disposable insufflation bulb and pilot balloon to provide a complete and economical single use system.
- **SAFETY:** Should local hospital policy require it, the battery pack is fully removable to permit safe, environmentally responsible disposal. Each carton is supplied with a battery storage tray to ensure safe separation of used battery cells.



| Description | Order Code |
|--|-------------------------------------|
| Rocket® Pro-Lux™ Sigmoidoscope Self-contained illumination handle, insufflation bulb, pilot balloon. Non-sterile, 25 units per carton | R53901 USL: 853901 |



North Island: 0800 658 814 South Island: 0800 804 546
customerservices@uslmedical.co.nz

nusurgix is a trading name of Rocket Medical plc