




With NRFit® connectors

- NRFit® SPINAL MANOMETERS**
 ensure connection to ISO80369-6 compliant spinal needles and syringes for the safe removal and measurement of cerebrospinal fluid (CSF).
- STACKABLE ASSEMBLY:** NRFit® Male/female tube fittings permit additional sections to be stacked, allowing pressures $>40\text{cmH}_2\text{O}$ to be measured when required.
- CLEAR, STURDY PLASTIC CONSTRUCTION** for reliable, secure assembly and ease of use.
- SUPPLIED COMPLETE** with 3-way stopcock for attachment to all NRFit® ISO80369-6 compliant spinal needles and syringes.
- READY FOR IMMEDIATE USE** for when time is precious.

Description	Order Code	NHSSC Code
Spinal Manometer - NRFit® Sterile, for single use, packed in cartons of 10 units	R55990-ISO-6	eDirect GBR000



Spinal Manometers have been used to provide safe, reliable, cost-effective CSF pressure measurement for over 30 years.

The unique 'stackable' design, original to Rocket Medical, has become a standard neurological diagnostic tool and can be found in emergency departments, neurological wards, medical wards and specialist neurosurgery centres throughout the world.

This original Rocket Medical design is now available with ISO80369-6 NRFit®-compliant connectors.

NRFit® is a trademark of GEDSA and is used with permission.

Rocket Medical Plc

Sedling Road, Washington, Tyne & Wear, NE38 9BZ

T: +44 (0) 191 419 6988 E: customerservices@rocketmedical.com W: www.rocketmedical.com

 SALES 384 Rev: 07
 Date: 2021-04-14
 Copyright© 2017-21