



NRFit®


CONTENTS:

- 1 x NRFit® 40cm stackable Spinal Manometer Set
- 1 x NRFit® 22G x 9cm QUINKE spinal needle
- 1 x 10ml NRFit® syringe
- 5 x 10cm x 10cm gauze swabs
- 1 x 60ml gallipot
- 1 x 8" skin prep sponge stick
- 1 x adhesive fenestrated drape, 75 x 90cm
- 1 x square drape 1m²
- 3 x 30ml Universal Containers with labels
- 1 x 20ml luer lock syringe
- 1 x microlance 21G /1.5" green needle (luer)
- 1 x microlance 25G/5/8" yellow needle (luer)

- **KIT INCLUDES Pajunk NRFit® SPINAL NEEDLE:** 22G x 9cm QUINKE spinal needle.
- **REDUCES RISK:** Combining all the necessary NRFit®-compatible components into a single kit significantly reduces the risk of mismatch or accidental connection to non-approved systems.
- **STURDY CLEAR CONSTRUCTION** for reliable, secure assembly and ease of use. Ready for immediate use for when time is precious.
- **TWO PART MANOMETER ASSEMBLY** allows CSF pressure monitoring between 3 - 40cmH₂O. Supplied with 3 way stopcock attached for connection to the NRFit® compliant spinal needle and syringe included in the kit.

Ordering Information	Code
Lumbar Puncture Tray with NRFit® Spinal Manometer Complete ward or department procedure kit for single patient use Packed in cartons 5 units	R55995-NRFIT

NRFit® is a trademark of GEDSA and is used with permission.