



NRFit®


CONTENTS:

- 1 x NRFit® 40cm H₂O stackable spinal manometer
- 1 x 10ml NRFit® slip syringe
- 5 x 10cm x 10cm gauze swabs
- 1 x 60mL pot
- 1 x Skin prep sponge stick
- 1 x 2-piece fenestrated drape, 75 x 90cm (75 x 45cm x 2)
- 1 x square drape 60 x 60 cm
- 3 x 30mL Universal Containers with labels
- 1 x 20mL luer lock syringe
- 1 x Microlance 21G /1.5" green needle (luer)
- 1 x Surguard 2 25G /38mm orange needle (luer)

- **CONVENIENT:** Combines the necessary* components for lumbar puncture into a single kit so that they are at hand when you need them.
*As needle choice is an important consideration for clinicians, a spinal needle is not included in the kit; this provides users with the flexibility to use an NRFit® needle of their choice.
- **ENHANCED SAFETY:** The manometer features NRFit® connections. These connections feature on devices used for the administration of fluids to neuraxial systems or in regional anaesthesia. NRFit® connections are incompatible with luer devices, which prevents the accidental misconnection of devices used for the administration of fluids intended for other bodily systems.
- **STURDY, CLEAR CONSTRUCTION** for reliable, secure assembly and ease of use.
- **TWO-PART MANOMETER ASSEMBLY** allows CSF pressure monitoring between 3 - 40cm H₂O. Supplied with an attached 3-way stopcock, for connection to the NRFit® syringe included in the kit, and an NRFit® spinal needle of your choice.

Ordering Information	Code	USL Code
Lumbar Puncture Tray with NRFit® Spinal Manometer Complete procedure kit, for single patient use. Packed in cartons of 5 units.	R55995-NRFIT	855995-NRFIT



North Island: 0800 658 814 South Island: 0800 804 546

NRFit® is a trademark of GEDSA and is used with permission.

Rocket Medical Ltd.

PO Box 646, Orewa, 0946, New Zealand
 T: +64 (0) 21 198 8488 E: nz@rocketmedical.com W: nz.rocketmedical.com

Sales: 586 Rev: 02
 Date: 2022-02-10
 Copyright© 2021-22