

Ultra-Fine™ LLETZ Loops & Ball Electrodes



Ultra-Fine™ LLETZ Loops were the world's first single use loop electrodes designed specially for LLETZ.

Ultra-Fine™ LLETZ Loops have been in continuous production for more than 27 years. They continue our tradition of providing exceptional value, reliability and quality

- **OVERCOME THE PROBLEMS** of high cost, cross infection risks and health and safety hazards commonly associated with reprocessing reusable LLETZ electrodes.
- **Ultra-Fine™ LLETZ LOOPS** are designed to deliver higher cutting performance with less thermal damage to the sample.
- **2.4mm (3/32") x 10cm shaft** to fit all standard UK, US and European diathermy systems.

Ultra-Fine™ LLETZ Range	Size		Code
Ultra-Fine™ LLETZ Loop, for single use, packed sterile in cartons 25 units	Small	14mm	R57800
	Medium	18mm	R57801
	Large	22mm	R57802
Ultra-Fine™ Ball Electrode. for single use, packed sterile in cartons 25 units	Standard	5mm	R57804
	Small	3mm	R57807

