

## Water Trap Sets for Use with CRAFT™ Suction Pumps

- **COMPLETE SET:** Provides the complete bottle and tubing set in a sterile pack for single patient use.
- **RAPID DETECTION OF POTENTIAL CONTAMINATION:** The clear glass bottle enables the rapid detection of aspirate, should it enter the trap. Designed to withstand >450mmHg vacuum.
- **RELIABLE CONNECTION:** 2m patient tube set with integral luer lock connector to ensure secure, reliable connection to the aspiration needle set.
- **SUITABLE FOR ALL CRAFT™ PUMPS:** Compatible with all versions of CRAFT™ Suction and CRAFT™ DUO-VAC™ oocyte aspiration pumps.

Description	Order Code
<b>Single Use Water Trap Set:</b> Suitable for CRAFT™ Suction and CRAFT™ DUO-VAC™ Pumps Supplied sterile, for single use and in cartons of 10 units.	<b>R57685</b>



**The Water Trap Sets are designed to protect the CRAFT™ suction pumps from contamination.**

Should any fluid be aspirated towards the pump it is retained in the trap.

The set should be changed immediately should any liquid be collected, and between patients.

